

SECTION**9****PARTS LIST 1310 , 1320 & 1321F**

Description	115V	220V
Cord Set – European	NA	1800500
Cord Set – USA	1800516	101990
Door Gasket	600002	600002
Door Handle	300142	300142
CIRCUIT BREAKER	1100505	1100505
Fan 1321F	2600502	2600502
Fan Blade 1321F	104073	104073
Heating Element, 1310	9570637	9570642
Heating Element, 1320	9570640	9570644
Heating Element 1321F	9570746	9570802
Knob, Main Temperature	X1000771	X1000771
Knob, Over Temperature	4450506	4450506
Main Temperature Controller	100026	100026
Motor 1321F	4880527	4880528
On/Off (I/O) Switch	103351	103351
Over Temperature Thermostat	1750615	1750615
Pilot Light, green	200021	200021
Pilot Light, red	200020	200020
Shelf Clips	200116	200116
Shelf, 1310	5080538	5080538
Shelf, 1320, 1321F	5080539	5080539
Thermometer	200063	200063
Thermometer Clip	5080865	5080865

UNIT SPECIFICATIONS

Weight	Shipping	Net
1310	45 lbs.	38 lbs.
1320	67 lbs.	60 lbs.

Dimensions	Exterior WxDxH	Interior WxDxH
1310	18X16.5X21.25"	12 X 12 X 14"
1320	21 X 17.5 X 23.5"	17 X 12 X 17"

Capacity	Cubic Feet
1310	1
1320	2

Temperature	Range	Uniformity
1310	Amb+5° to 200°C	+4°C
1320	Amb+5° to 200°C	+4°C