

UNIT SPECIFICATIONS

These units are 120 volt or 230 volt. Please refer to the unit data plate for individual electrical specifications.

Table 8. Weight

Model	Shipping	Net
1525	120 lbs.	83 lbs.
1535	204 lbs.	162 lbs.
1555	300 lbs.	195 lbs.
1545	275 lbs.	158 lbs.
1565	550 lbs.	316 lbs.

Table 9. Dimensions

Model	Exterior WxDxH	Interior WxDxH
1525	21 x 21 x 25 in.	15 x 15 x 15 in.
1535	30 x 30 x 32 in.	24 x 24 x 20 in.
1555	42.25 x 20 x 37 in.	36 x 20 x 26 in.
1545	25 x 27 x 37 in.	19.25 x 20 x 26 in.
1565	25 x 27 x 72 in.	19.25 x 20 x 26 in.*

Note: *Each chamber.

Table 10. Capacity

Model	Cubic Feet
1525	2.0
1535	6.7
1555	11.3
1545	5.8
1565	5.8 (each chamber)

Note: *Each chamber.

Table 11. Temperature

Model	Range	Uniformity	Stability
1525	Amb. +5° to 70°C	+.25° @ 37°C	+.1°C
1535	Amb. +5° to 70°C	+.25° @ 37°C	+.1°C
1555	Amb. +5° to 70°C	+.25° @ 37°C	+.1°C
1545	Amb. +5° to 70°C	+.25° @ 37°C	+.1°C
1565	Amb. +5° to 70°C	+.25° @ 37°C	+.1°C

Table 12. Power

Model	Voltage	Amperage
1525	115 Volts	3.0 Amps
1525-2	220 Volts	1.5 Amps
1535	115 Volts	5.0 Amps
1535-2	220 Volts	2.5 Amps
1555	115 Volts	7.0 Amps
1555-2	220 Volts	4.0 Amps
1545	115 Volts	.5 Amps
1545-2	220 Volts	2.5 Amps
1565	115 Volts	.5 Amps
1565-2	220 Volts	2.5 Amps