



Overall dimensions (L x W x H):	17 x 11 x 7.75" (43.2 x 27.9 x 19.7cm)
Tray dimensions (L x W):	11 x 7.75" (27.9 x 19.7cm)
Electrical (50/60 Hz):	120 volts, 5 amps, 450 watts 230 volts, 5 amps, 450 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Temperature range:	ambient +5°C to 65°C
Temperature uniformity:	±0.5°C at 37°C
Speed range:	100 to 1200rpm
Speed accuracy:	±2% of set speed up to 999rpm ±5% 1000 to 1200rpm
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	0.125" (3mm)
Capacity:	4 microplates or 2 Micro-Tube Racks
Controls:	see page 4
Tray material:	aluminum
Ship weight:	30lbs (13.6kg)