

MINI SHAKER SPECIFICATIONS



Overall dimensions (L x W x H):	17 x 11 x 4" (43.2 x 27.9 x 10.2cm)
Tray dimensions (L x W):	11.75 x 8.75" (29.9 x 22.2cm)
Electrical (50/60 Hz):	120 volts, 5 amps, 25 watts 230 volts, 5 amps, 25 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Speed range:	100 to 1200rpm
Speed accuracy:	±2%
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	0.125" (3mm)
Maximum weight capacity:	~ 8lbs (3.6kg), up to 1000rpm ~ 5lbs (2.3kg), over 1000rpm
Controls:	see page 4
Tray material:	aluminum
Ship weight:	25lbs (11.4kg)

Mini Shaker with plasticware