



Overall dimensions (L x W x H):	17 x 11 x 5" (43.2 x 27.9 x 12.7cm)
Platform dimensions (L x W):	12.75 x 10" (32.4 x 25.4cm)
Electrical (50/60 Hz):	120 volts AC, 5 amps, 25 watts 230 volts AC, 5 amps, 25 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Speed range:	1 to 50rpm*
Speed accuracy:	±1rpm
Tilt angle:	0 to 15°*
Weight capacity:	10lbs (4.5kg)**
Timer:	digital, 1 second to 9999 minutes (increased in 1 second increments)
Controls:	see page 4
Ship weight:	15.5lbs (7kg)

* Maximum speed/tilt angle may vary with heavy or unbalanced loads.

** Centered on tray.

VWR Rocker with glassware and plasticware accessories