

SPECIFICATIONS - ANALOG MULTI-TUBE VORTEXER

Overall dimensions (L x W x H): 9.5 x 15.1 x 16" (24.1 x 38.4 x 40.6cm)

Electrical (50/60 Hz):

120V - 0.8 amps, 100 watts

230V - 0.4 amps, 100 watts

Fuses: 120V - 5mm x 20mm, 5 amp, 250V quick acting

230V - 5mm x 20mm, 5 amp, 250V quick acting

Speed range: 1200 to 2400rpm

Accuracy: +/- 25rpm

Orbit: 3.6mm

Weight capacity: 10lbs (4.5kg)

Timer: mechanical timer runs 0 to 60 seconds

Controls: rocker switch: run, standby, timer, speed knob,

variable 1 to 10 dial markings, time knob, vari-

able 1 to 10 dial markings

Ship weight: 45lbs (20.4kg)

