



## BenchRocker™ 3D Variable Speed Gyrotory Rocker (3D)



WITH  
**Insta  
TILT™**

- + Adjustable speed & tilt for a wide range of applications
- + Large 14x12" work surface
- + Insta-TILT™ technology for instant adjustment of tilt angle
- + Cold room and incubator safe



Optional Stacking Platform (clearance: 3")

CE

The BenchRocker 3D is designed to provide the precise speeds and tilt angles required for a broad range of molecular and biological mixing applications. The user can adjust both parameters to match the vessel size and the volume of liquid being mixed, yielding optimum results. Proprietary InstaTilt™ technology means no tools, no disassembly. ***Simply grip the platform with two hands and move to the desired position.***

This model includes both flat and dimpled (autoclavable) mats for mixing a wide variety of tubes and other vessels. Stacking platforms can be purchased as accessories, effectively doubling the usable surface area of each model.

BenchRocker 3D features a maintenance free, brushless motor, capable of loads up to 4.4 lbs and is safe for use in cold rooms and incubators.

### Ordering Information:

| Item No.       | Description:   |
|----------------|--|
| B3D2300*       | BenchRocker 3D with flat and dimpled mats (12x12")                 |
| BR1000-STACK   | Stacking platform with flat mat (12x12")                           |
| BR1000-STACK-D | Stacking platform with dimpled mat (12x12")                        |
| B0718          | Bungee cord, used for securing vessels, pack of 4                  |
| BR2000-SP      | STACK-PLUS Extender set, adds 1.25" clearance to stacking platform |

\*115V with US plug. For 230v (EU plug), please add (-E).

### Technical Data:

|                        |                               |
|------------------------|-------------------------------|
| Speed Range:           | Variable, 2 to 30 rpm         |
| Tilt Angle:            | Variable, 0 to $\pm 30^\circ$ |
| Platform Size:         | 12 x 12 in. / 30 x 30 cm      |
| Load Capacity:         | 5 lbs (2.3kg)                 |
| Operating Temp. Range: | +4°C to +65°C                 |
| Dimensions:            | 12 x 12 x 8 in.               |
| (W x D x H)            | 30 x 30 x 20 cm               |
| Weight:                | 7.5 lbs (3.4kg)               |
| Electrical:            | 115 or 230V, 50-60Hz, 50W     |
| Warranty:              | 2 Years                       |