



Fisher Scientific™ Analog Vortex Mixer

Variable speed, analog control that allows low rpm start-up for gentle shaking or high-speed mixing for vigorous vortexing of samples

Catalog No. 02-215-414

Includes: Includes both cup head and 3" head with cover. 72 inch (183cm), 3-wire cord and plug.



Specifications

Speed	300 to 3200rpm
Height (English) Exterior	5 in.
Length (English) Exterior	7.97 in.
Width (English) Exterior	4.75 in.
Height (Metric) Exterior	12.7cm
Length (Metric) Exterior	20.2cm
Width (Metric) Exterior	12.1cm
Shipping Weight (English)	10 lbs
Shipping Weight (Metric)	4.5kg
Voltage	120V
Amperage	1.2 amps
Plug Type	US
Warranty and Service Offering	5 year limited warranty on parts and labor.
Includes	Includes both cup head and 3 inch head with cover. Supplied with a 72 inch, 3-wire cord and plug.
Orbit	4.9mm
Wattage	120V
Weight (English)	10 lb.
Weight (Metric)	4.5kg
Dimensions (L x W x H) Exterior	7.97 x 4.75 x 5 in. (20.2 x 12.1 x 12.7 cm)



A Thermo Fisher Scientific Brand

Electrical Requirements	120V, 1.2 amps, 150 watts
Length (Metric)	20.2cm
Length (English)	7.97 in.
Width (English)	4.75 in.
Height (English)	5 in.
Width (Metric)	12.1cm
Height (Metric)	12.7cm
Warranty	5 year limited warranty on parts and labor.