


## 211DS Digital Shaking Incubator



S2031-13 with K566



- Integrated digital orbital shaker
- Exclusive SmartCheck temperature control system
- Broad temperature ranges to meet a wide variety of applications
- Stacking design saves space

### Ordering Information

Cat#	Description
I5211-DS	<b>211DS</b> Shaking Incubator, 120V
I5211-SA	Stacking adapter, required when stacking incubators
I5230	Flask clamp platform, pre drilled (see clamps below)
I5231	Flat platform with non-slip rubber mat (30 x 30cm)
I5231-D	Double flat platform with non-slip rubber mat
S2031-13	Universal spring loaded clamp platform for flasks, bottles, etc.
S2040-85	Flask clamp, 125mL (max 16)
S2040-99	Flask clamp, 250mL (max 9)
S2040-09	Flask clamp, 500mL (max 5)
S2040-01	Flask clamp, 1000mL (max 4)

\*To order 230V units add -230V to the end of the catalog number. 230V unit includes -EU and -UK cords  
\*\*Other size flask clamps are available. Call for details.

### Specifications

Temp Range	Ambient +5° to 80°C
Temp Uniformity	±0.25°C
Temp Accuracy	±0.1°C
Shaker speed	20 - 300 rpm
Timer	0-99 min or continuous
Shaker orbit	19 mm
Exterior dims	42.5 x 55 x 58 cm/ 16.75 x 21.75 x 23 in
Chamber dims	34.3 x 37.5 x 38.1 cm/ 13.5 x 14.25 x 15 in
Chamber volume	1.7 cu ft/49 liters
Weight	83 lb/37.7 kg



## 311DS Digital Shaking Incubator

- Exclusive SmartCheck temperature control system
- Broad temperature range to meet a wide variety of applications
- Gasketed door with large, integral viewing area

### Ordering Information

Cat#	Description
I-5311-DS*	<b>311DS</b> Shaking Incubator, 120V
I5322	Additional shelf
I5331	Flat platform with non-slip rubber mat (39 x 33 cm)
I5330	Shaking platform pre-drilled for flask clamps (flask clamps sold separately, see 211DS)
I5330-2000	Platform with 4 x 2L flask clamps
I5330-1000	Platform with 6 x 1L flask clamps
I5330-500	Platform with 8 x 500mL flask clamps
I5330-250	Platform with 12 x 250mL flask clamps
I5330-125	Platform with 20 x 125mL flask clamps
I5330-50	Platform with 30 x 50mL flask clamps

\*To order 230V units add -230V to the end of the catalog number. 230V units include -EU and -UK cords.  
\*\*Rack holders to fit onto universal platform are available. Call for details.



### Specifications

Temp Range	Ambient +5° to 80°C
Temp Uniformity	±0.25°C
Temp Accuracy	±0.1°C
Shaker speed	20 - 300 rpm
Timer	0-99 min or continuous
Shaker orbit	19 mm
Exterior dims	57.4 x 54.4 x 63.5 cm/ 22.6 x 21.4 x 25 in
Chamber dims	44.2 x 39.6 x 40.6 cm/ 17.4 x 15.6 x 16 in
Chamber volume	2.5 cu ft/71 liters
Weight	140 lb/63.6 kg

