

# Specifications

---

## Trays

---

### Tray Dimensions

Large - 15 1/2" x 9 3/16"  
(39.4 x 23.3 cm)

Small - 15 1/2" x 6"  
(39.4 x 15.2 cm)

---

### Number of Trays

The sterilizer is supplied with two trays, 1 large and 1 small, and one rack.

---

## Dimensions

---

### Capacity

34 Liters (Vessel) (1.2 cubic feet)  
7.0 Liters (Reservoir)

---

### Overall Dimensions

W x D x H Inches(cm)  
22 7/16" x 24 1/2" x 17"  
(57.0 x 62.2 x 43.2)  
Chamber (Dia. x Depth)  
12" X 18"  
(30.5 x 45.7)

---

## Shipping Weight

85 lbs. (38.5 kg)

---

## Electrical

---

### Heater Power

1500 W

---

### Power requirements

120V, 50/60 Hz, 12.5 Amps  
220V~, 50/60 Hz, 6.5 Amps  
100V, 50/60 Hz, 15 Amps

---

## Standards

ASME  
UL/cUL to UL3101-1

---

---

## Reservoir Water

Distilled or low quality  
deionized  
Max. 1.5 Megohm/cm  
Min. .5 Megohm/cm

---

## Mounting

Benchtop only

---

## Controls

### Fixed Cycles

- Liquid
- Packs
- Wrapped
- Unwrapped

### Optional Cycle



### Note

You can examine the pressure vessel rating plate for the sterilizer chamber by opening the panel on the back cover of the sterilizer, then removing the small plug of insulation on the back of the sterilizer chamber.